

Arcadia14 Medical Supplies Ltd.  
Hill House, Office 117, 210 Upper Richmond  
Road, Putney, London SW15 6NP

**Phone:** +44 203 488 3521

**Enquires:** [info@arcadia14.co.uk](mailto:info@arcadia14.co.uk)


**Sales:** [sales@arcadia14.co.uk](mailto:sales@arcadia14.co.uk)



**Dental**  
Instruments Catalogue

[www.arcadia14.co.uk](http://www.arcadia14.co.uk)


# Table Of Contents




Extracting Forceps 01-31



Root Elevators 32-51




Filling Instruments 52-56



Dental Probes 57-60



Periodontia Instruments & Mallets 61-68



Curettes & Scalers 69-75

# Table Of Contents


Modelling & Measuring  
Instruments 76-78




Wax / Plaster Knives 79-80



Diagnostic (Explorers) &  
Endodontic Instruments 81-86



Excavators 87-88



Handles For Mouth Mirrors  
& Mouth Mirrors 89-90



Spatulas 91-94



# Table Of Contents



Needle Holders &  
Saliva Ejectors

95-100



Syringes

101-104



Dental Tweezers

105-109



Rubber Dam Clamps

110-112



Matrix Retainers &  
Matrix Bands

113-114



Crown Removers, Scissors,  
Forceps, Amalgum Guns & Carriers

115-117

# Table Of Contents

Copper Ring Removing Pliers & Paper Articulating Forceps 118




Rongeurs 119-122



Bone Elevators, Files & Chisels 123-124



Pliers For Orthodontics & Prosthetics 125-135



Impression Trays 136-139



Articulators 140-141



Hollowares 142-144





## SECTION **01**

**Extracting Forceps - English Pattern**

**Page 01-31**

**[www.arcadia14.co.uk](http://www.arcadia14.co.uk)**

# Extracting Forceps - English Pattern



**1001**  
Fig 8  
lower premolars



**1002**  
Fig 13  
lower premolars



**1003**  
Fig 17  
upper molars, right



**1004**  
Fig 18  
upper molars, left




**1005**  
Fig 18 A  
upper molars





**1006**  
Fig 19  
upper wisdoms



**1007**  
Fig 20  
lower wisdoms



**1008**  
Fig 21  
lower molars



**1009**  
Fig 22  
lower molars



**1010**  
Fig 22 A  
Albrecht  
lower molars

# Extracting Forceps - English Pattern



**1011**  
Fig 33 A  
lower roots



**1012**  
Fig 33 B  
lower roots



**1013**  
Fig 33 L  
lower roots



**1014**  
Fig 33 L  
lower roots




**1015**  
Fig 33 S  
lower roots, for children



**1016**  
Fig 44  
upper roots




**1017**  
Fig 45  
lower roots



**1018**  
Fig 12  
upper roots



**1019**  
Fig 51 A  
upper roots



**1020**  
Fig 51 L  
upper roots

**"Relax" - Extracting Forceps**



**Relax  
1021**  
Fig 51  
upper roots



**Relax  
1022**  
Fig 51 A  
upper roots



**Relax  
1023**  
Fig 51 L  
upper roots



**Relax  
1024**  
Fig 67  
upper wisdoms



**Relax  
1025**  
Fig 67 A  
upper wisdoms



**1026**  
Fig 37  
Upper incisors



**1027**  
Fig 39  
Upper molars



**1028**  
Fig 39 R  
Upper molars, right



**1029**  
Fig 39 L  
Upper molars, left



**1030**  
Fig 515  
Upper roots

# Extracting Forceps For Children, English Pattern




**Klein  
1031  
Fig 6  
lower molars**



**Klein  
1032  
Fig 137  
upper incisors**



**Klein  
1033  
Fig 5  
lower incisors**



**Klein  
1034  
Fig 139  
uper premolars**




**1035**  
Fig 43 A  
Fractured or decayed  
upper teeth



**1036**  
Fig 44  
Upper roots



**1037**  
Fig 45  
Lower roots



**1038**  
Fig 46  
Lower roots



**1039**  
Fig 46 L  
Lower roots

# Extracting Forceps - English Pattern



**1040**  
Fig 48  
Upper roots



**1041**  
Fig 51  
Upper roots



**1042**  
Fig 49  
Upper roots



**1043**  
Fig 51 L  
Upper roots



**1044**  
Fig 51 A  
Upper roots



**1045**  
Fig 51 LX  
Upper roots



**1046**  
Fig 52  
Upper roots and  
canines



**1047**  
Fig 54  
For separating  
upper molars



**1048**  
Fig 55  
For separating  
upper molars



**1049**  
Fig 56  
For separating  
lower molars

# Extracting Forceps - English Pattern





**1055**  
Fig 66 L  
Upper molars, left



WITH SERRATED  
TIPS

**1056**  
Fig 66 R  
Upper molars, right



WITH SERRATED  
TIPS

**1057**  
Fig 67  
Upper third molars



WITH SERRATED  
TIPS

**1058**  
Fig 67 A  
Upper third molars



WITH SERRATED  
TIPS

**1059**  
Fig 68  
Lower roots

# Extracting Forceps - English Pattern



**1060**  
Fig 67½ L  
Upper third molars, left



**1061**  
Fig 67½ R  
Upper third molars, right



**1062**  
Fig 73  
Lower molars



**1063**  
Fig 74  
Lower roots



**1064**  
Fig 74 M  
Lower roots



**1065**  
Fig 74 N  
Lower roots

WITH SERRATED  
TIPS



**1066**  
Fig 75  
Lower canines and  
premolars

WITH SERRATED  
TIPS



**1067**  
Fig 79  
Lower third molars

WITH SERRATED  
TIPS



**1068**  
Fig 79 A  
Lower third molars

WITH SERRATED  
TIPS



**1069**  
Fig 79 C  
Lower third molars

WITH SERRATED  
TIPS

# Extracting Forceps - English Pattern



**1070**  
Fig 86 W  
Lower third molars



**1071**  
Fig 86 1/2 L  
Lower third molars, left



**1072**  
Fig 86 1/2 R  
Lower third molars, right



**1073**  
Fig 86 C  
Lower molars



**1074**  
Fig 87  
Lower decayed molars  
or broken-down crowns



WITH SERRATED  
TIPS

**1075**  
Fig 89  
Upper molars, right



WITH SERRATED  
TIPS

**1076**  
Fig 90  
Upper molars, left



WITH SERRATED  
TIPS

**1077**  
Fig 99  
Broken-down lower molars



WITH SERRATED  
TIPS

**1078**  
Fig 107  
Upper canines

# Extracting Forceps - English Pattern



**1079**  
Fig 121  
Upper third molars



**1080**  
Fig 145  
Lower roots



**1081**  
Fig 151½  
Upper roots



**1082**  
Fig 166  
Universal extracting forceps for  
upper teeth, in particular molars  
and premolars



**1083**  
Fig 167  
Universal extracting forceps for  
lower teeth, in particular molars  
and premolars



**1084**  
Fig 168  
Upper canines



**1085**  
Fig 169  
Lower canines



**1086**  
Fig 346  
Lower roots



**1087**  
Fig 349  
Upper roots



**1088**  
Fig 359  
Lower roots

# Extracting Forceps - Mead Pattern



**1089**  
Fig 1  
Universal extracting forcep for  
upper incisors and canines



**1090**  
Fig 2  
Upper molars



**1091**  
Fig 3  
Lower incisors, premolar  
and roots



**1092**  
Fig 4  
Lower molars



**1093**  
Fig 1  
Upper incisors  
and canines



**1094**  
Fig 2  
Upper incisors  
and canines



**1095**  
Fig 4  
Lower incisors  
and canines



**1096**  
Fig 7  
Upper premolars



**1097**  
Fig 13  
Lower premolars

## Extracting Forceps With Anatomically Shyped Handle



**1098**  
Fig 17  
Upper molars, right



WITH SERRATED  
TIPS

**1099**  
Fig 18  
Upper molars, left



WITH SERRATED  
TIPS

**1100**  
Fig 21  
Lower molars



WITH SERRATED  
TIPS

**1101**  
Fig 22  
Lower molars



WITH SERRATED  
TIPS

**1102**  
Fig 33  
Lower roots



**1103**  
Fig 45  
Lower roots



WITH SERRATED  
TIPS



WITH SERRATED  
TIPS

**1104**  
Fig 51 A  
Upper roots



WITH SERRATED  
TIPS

**1105**  
Fig 67 A  
Upper third  
molars



WITH SERRATED  
TIPS

**1106**  
Fig 74  
Lower roots



WITH SERRATED  
TIPS

**1107**  
Fig 79  
Lower third  
molars



WITH SERRATED  
TIPS

**1108**  
Fig 86 C  
Lower molars



WITH SERRATED  
TIPS

**1109**  
Fig 359  
Lower roots

# Extracting Forceps With Anatomically Shyped Handle



**1110**  
Fig 89  
Upper molars, right



**1111**  
Fig 90  
Upper molars, left



**1112**  
Fig 145  
Lower roots

www.arcadia14.co.uk



**1113**  
Fig 151 1/2  
Upper roots



**1114**  
Fig 349  
Upper roots



**1115**  
Fig 5  
Lower third molars



**1116**  
Fig 6  
Separating forcep



**1117**  
Fig 105  
Upper molars



**1118**  
Fig 16  
Lower molars



**1119**  
Fig 17  
Lower molars

# Extracting Forceps - American Pattern



**1120**  
Fig 18 L  
Upper molars, left



**1121**  
Fig 18 R  
Upper molars, right



**1122**  
Fig 23  
Lower molars



**1123**  
Fig 53 R  
Upper molars, right



**1124**  
Fig 53 L  
Upper molars, left



**1125**  
Fig 65  
Upper incisors and roots



**1126**  
Fig 69  
Upper and lower roots



**1127**  
Fig 88 R  
Upper molars, right



**1128**  
Fig 88 L  
Upper molars, left



**1129**  
Fig 101  
Deciduous teeth

## Extracting Forceps - American Pattern



**1130**  
Fig 150  
Upper incisors, premolars, roots



**1131**  
Fig 150 S  
Upper teeth (for children)



**1132**  
Fig 151  
Upper incisors, premolars, roots



**1133**  
Fig 151S  
Lower teeth (for children)



**1134**  
Fig 222  
Lower third molars




**1135**  
Fig 1  
upper incisors and bicusps



**1136**  
Fig 1 B  
upper roots, anterior



**Physick**  
**1137**  
Fig 5  
lower wisdoms



**1138**  
Fig 6  
lower molars, splitting forceps



**1139**  
Fig 10 H  
upper molars

# Extracting Forceps For Children - English Pattern



**1140**  
Fig 135  
Lower premolars

WITH SERRATED  
TIPS



**1141**  
Fig 225  
Lower molars



**1142**  
Fig 305  
Upper incisors  
and roots



**1143**  
Fig 33 S  
Lower roots

WITH SERRATED  
TIPS



**1144**  
Fig 38  
Lower incisors  
and canines

WITH SERRATED  
TIPS



**1145**  
Fig 1  
Upper incisors



**1146**  
Fig 2  
Upper premolars



**1147**  
Fig 3  
Upper molars



**1148**  
Fig 4  
Upper roots



**1149**  
Fig 5  
Lower incisors



**1150**  
Fig 6  
Lower molars



**1151**  
Fig 7  
Lower roots

# Root Splinter Extracting Forceps



**1152**  
Fig 1  
Lower roots



**1153**  
Fig 2  
Upper roots



**1155**  
Fig 2  
Roots fragments



**1154**  
Fig 1  
Roots fragments

## SECTION **02**


ROOT ELEVATORS, FILLING  
INSTRUMENTS & DENTAL PROBES

Page 32-60

[www.arcadia14.co.uk](http://www.arcadia14.co.uk)

lower jaw, left side

lower jaw, right side



**1200**



**1201**

Root elevator, spiral-shaped, for extremely facilitated extraction of the molars in the lower jaw. The slender, conical points allow deep penetration between tooth and gum. Owing to the hook-shaped working parts, slipping of the elevator is avoided.

## 02 SECTION

### Root Elevators



Bein



**Bein  
1202  
Fig 1**



**Bein  
1203  
Fig 2**



**Bein  
1204  
Fig 0**



**Bein  
1205  
Fig 3**



**Bein  
1206  
Fig 4**



**Bein  
1207  
Fig 5**



**Bein  
1208  
Fig 6**



**Seldin**



**Seldin  
1209  
Fig 1L**



**Seldin  
1210  
Fig 1R**



**Seldin  
1211  
Fig 2**



**Seldin  
1212  
Fig 3**



**Seldin  
1213  
Fig 4L**



**Seldin  
1214  
Fig 4R**

## 02 SECTION

### Root Elevators



**Cryer  
1215  
Fig 3**



**Cryer  
1216  
Fig 4**



**Cryer  
1217  
Fig 12**



**Cryer  
1218  
Fig 27**



**Cryer  
1219  
Fig 28**



**Cryer  
1220  
Fig 27 S**



**Cryer  
1221  
Fig 28 S**



**Apical**



**Apical  
1222  
Fig 1**



**Apical  
1223  
Fig 2**




**Apical  
1224  
Fig 3**



**Apical  
1225  
Fig 4**

## 02 SECTION

### Root Elevators





**1231**  
 Fig 1



**1232**  
 Fig 2



**1233**  
 Fig 3



**Flohr  
 1234**  
 Fig 1



**Flohr  
 1235**  
 Fig 2



**Flohr  
 1236**  
 Fig 3

## 02 SECTION

### Root Elevators



**Kopp  
1237  
Fig 1**



**Kopp  
1238  
Fig 2**



**Kopp  
1243  
Fig 1**



**Kopp  
1244  
Fig 2**



**Kopp**



**Kopp  
1239  
Fig 3**



**Kopp  
1240  
Fig 4**



**Kopp**



**Kopp  
1241  
Fig 5**



**Kopp  
1242  
Fig 6**



**Kopp  
1245  
Fig 3**



**Kopp  
1246  
Fig 4**



## 02 SECTION

### Root Elevators



**Warwick-James  
1253  
Fig 1**



**Warwick-James  
1254  
Fig 2**



**Warwick-James  
1255  
Fig 3**



**London Hospital  
1256  
Fig 1**




**London Hospital  
1257  
Fig 2**



**London Hospital  
1258  
Fig 3**




**Ott  
1259**



## 02 SECTION

### Root Elevators



**Mead  
1264  
MD 1**



**Mead  
1265  
MD-5 A**



**Mead  
1266  
MD-6 A**



**Mead  
1267  
MD-7 A**



**Mead  
1268  
MD-8**



**Lecluse  
1269**



**Chompret  
 1270**  
 Fig 1



**Chompret  
 1271**  
 Fig 2



**Chompret  
 1272**  
 Fig 3



**Chompret  
 1273**  
 Fig 4



**Apical  
 1274**



**1275**  
Fig 1



**1276**  
Fig 2



**1277**  
Fig 3



**Heidbrink**  
**1278**  
Fig 1



**Heidbrink**  
**1279**  
Fig 2



**Heidbrink**  
**1280**  
Fig 3



**Winter  
1281  
Fig 1L**




**Winter  
1282  
Fig 1R**



**Winter  
1283  
Fig 2L**



**Winter  
1284  
Fig 2R**



**Winter**



**Winter  
1285  
Fig 11 L**



**Winter  
1286  
Fig 11 R**



**Winter  
1287  
Fig 12 L**



**Winter  
1288  
Fig 12 R**



**Winter  
1289  
Fig 14 L**



**Winter  
1290  
Fig 14 R**

## 02 SECTION

### Root Elevators



**Barry  
1291  
Fig 1**



**Barry  
1292  
Fig 2**



**Barry  
1293  
Fig 3**



**Barry  
1294  
Fig 4**



**Barry**



**Barry  
1295  
Fig 11**



**Barry  
1296  
Fig 12**



**Barry  
1297  
Fig 17**



**Barry  
1298  
Fig 18**



**Cryer  
1300  
Fig 26**



**Seldin  
1301  
Fig 27**



**Cryer  
1302  
Fig 28**



**Cryer  
1303  
Fig 29**



**Heidbrink  
1304  
Fig 22**



**Heidbrink  
1305  
Fig 23**



**1299**

## 02 SECTION

### Root Elevators



**Kopp  
1307  
Fig 31**



**Kopp  
1308  
Fig 32**



**Kopp  
1309  
Fig 33**



**Kopp  
1310  
Fig 34**



**Kopp  
1311  
Fig 35**



**Kopp  
1312  
Fig 44**



**1306**



**Heidbrink**  
**1314**  
Fig 10



**Winter**  
**1315**  
Fig 11



**Winter**  
**1316**  
Fig 12



**Barry**  
**1317**  
Fig 17



**Barry**  
**1318**  
Fig 18



**1319**  
Fig 20



**1313**

## 02 SECTION

### Root Elevators

Root Elevators according to Dr. Sader,  
with interchangeable working tips




1321





1320



1322










**1443**  
Heidemann  
Fig 1



**1444**  
Heidemann  
Fig 2



**1445**  
Heidemann  
Fig 3



**1446**



**1447**  
Fig 146




**1448**  
Fig 147



**1449**  
Fig 148




**1450**  
Fig 149



**1451**  
Fig 158




**1452**  
Fig 159



**1453**  
Woodson  
Fig 160

# Periodontal Curettes and Filling Instruments



## 02 SECTION

### Dental Probes

With Point Set In



1470



1471-1



1472-2



1473-3



1474-4



1475-5



1476-6



1477-7



1478-8



1479-9



1480-10




1481-11



1482-2  
Woodson

## Burnishers and Pocket Probes Scalers

### Polishing Instruments





- New marking process.
- Laser system marking.
- No more paint peeling off.
- Permanent marking.
- Antiglare, low depth marking.
- Maximum tip resistance against breaking.



**1497**  
Fig CP-C




**1498**  
Fig CP-11



**1499**  
Fig CP-12



**1500**  
Fig CP-12S  
O.M.S



## SECTION **03**

PERIODONTIA INSTRUMENTS,  
MALLETS, CURETTES & SCALERS

Page 61-75

[www.arcadia14.co.uk](http://www.arcadia14.co.uk)

**Files**

**Chisels**



**1507**  
Sugarman  
Fig. 1S/2S



**1508**  
Sugarman  
Fig. 3S/4S



**1509**  
Schlugar  
Fig. 9/10



**1510**  
Buck  
Fig. 11/12



**1511**  
Ochsenbein  
Fig. 3




**1512**  
Ochsenbein  
Fig. 4



**1513**  
Rhodes  
Back - Action  
Fig. 36/37

### Files



**Gingivectomy knives**



**1521**  
Goldman-Fox  
Fig. 7



**1522**  
Goldman-Fox  
Fig. 8



**1523**  
Goldman-Fox  
Fig. 9



**1524**  
Goldman-Fox  
Fig. 11



**1525**  
Buck  
Fig. 3/4



**1526**  
Buck  
Fig. 5/6

### Curettes and Scalers



**1527**  
Curette  
Fig. HU 1/2



**1528**  
Columbia  
Fig. 2R/2L



**1529**  
Columbia  
Fig. 4R/4L



**1530**  
Columbia  
Fig. 13/14



**1531**  
Fig. 411/412




**1532**  
Fig. 413/414



**1533**  
Zerfing-chisel

**Gingivectomy knives**



### Curettes and Scalers



**1540**  
Curette  
Fig. HU 1/2



**1541**  
Columbia  
Fig. 2R/2L



**1542**  
Columbia  
Fig. 4R/4L



**1543**  
Columbia  
Fig. 13/14



**1544**  
Fig. 413/412




**1545**  
Fig. 411/414




**1546**  
Zerfing-chisel


**Bone Chisels and Gouges**




### Bone Chisels and Gouges



### Martelli chirurgici Mallet





**1576**  
Gracey  
Fig. 1/2



**1577**  
Gracey  
Fig. 3/4



**1578**  
Gracey  
Fig. 5/6



**1579**  
Gracey  
Fig. 7/8



**1580**  
Gracey  
Fig. 9/10



**1581**  
Gracey  
Fig. 11/12




**1582**  
Gracey  
Fig. 13/14



**1583**  
Gracey  
Mesial Distal  
Fig. 11/14



**1584**  
Gracey  
Mesial Distal  
Fig. 12/13



**1585**  
Gracey  
AF/type  
Fig. 1/2



**1586**  
Gracey  
AF/type  
Fig. 3/4



**1587**  
Gracey  
AF/type  
Fig. 5/6



**1588**  
Gracey  
AF/type  
Fig. 7/8



**1589**  
Gracey  
AF/type  
Fig. 11/12



**1590**  
Gracey  
AF/type  
Fig. 13/14



**1591**  
Gracey  
AF/type  
Fig. 15/16

# Gracey Curettes, AF, MF Type



**BH-1592**  
Gracey  
AF/type



**BH-1593**  
Gracey  
AF/type  
Fig. 1/2



**BH-1594**  
Gracey  
AF/type  
Fig. 3/4



**BH-1595**  
Gracey  
AF/type  
Fig. 5/6



**BH-1596**  
Gracey  
AF/type  
Fig. 7/8



**BH-1597**  
Gracey  
AF/type  
Fig. 11/12



**BH-1598**  
Gracey  
AF/type  
Fig. 13/14



**BH-1599**  
Gracey  
AF/type  
Fig. 15/16



**BH-1600**  
Gracey  
MF/type  
Fig. 1/2



**BH-1601**  
Gracey  
MF/type  
Fig. 3/4



**BH-1602**  
Gracey  
MF/type  
Fig. 5/6



**BH-1603**  
Gracey  
MF/type  
Fig. 7/8




**BH-1604**  
Gracey  
MF/type  
Fig. 11/12



**BH-1605**  
Gracey  
MF/type  
Fig. 13/14



**BH-1606**  
Gracey  
MF/type  
Fig. 15/16



**1607**  
Gracey  
MF/type  
Fig. 1/2



**1608**  
Gracey  
MF/type  
Fig. 3/4



**1609**  
Gracey  
MF/type  
Fig. 5/6



**1610**  
Gracey  
MF/type  
Fig. 7/8



**1611**  
Gracey  
MF/type  
Fig. 11/12



**1612**  
Gracey  
MF/type  
Fig. 13/14




**1613**  
Gracey  
MF/type  
Fig. 15/16



**1614**  
Goldman - Fax  
Fig. 1



**1615**  
Goldman - Fax  
Fig. 2




**1617**  
Goldman - Fax  
Fig. 3



**1618**  
Goldman - Fax  
Fig. 4



**1619**  
Goldman - Fax  
Fig. 21




**1620**  
Goldman - Fax  
Fig. 5



**1621**  
Mc Call  
Fig. 13S/14S



**1622**  
Younger-Good  
Fig. 7/8





**1635**  
Hemingway  
Fig. 1



**1636**  
Hemingway  
Fig. 2



**1637**  
Hemingway  
Fig. 3



**1638**  
Lucas  
Fig. 85



**1639**  
Lucas  
Fig. 86





**1640**  
Lucas  
Fig. 87

## SECTION **04**

MODELLING INSTRUMENTS, WAX KNIVES/PLASTER  
DIAGNOSTIC, EXPLORER, ENDODONTIC  
INSTRUMENTS & EXCAVATORS

Page 76-88

[www.arcadia14.co.uk](http://www.arcadia14.co.uk)





**1662**  
Boley




**1663**  
Zuricher  
130 mm



**1664**  
Vernier  
12 cm



**1665**  
Vernier




**1666**  
Carol



**1667**  
Castroviejo's  
Scleral Marker



**1668**  
Iwanson



**1669**  
Albee



**1700**  
18 cm  
7 1/8"



**1701**  
18 cm  
7 1/8"



**1702**  
17 cm  
6 3/4"  
Gritman



**1703**  
17 cm  
6 3/4"  
Medicon



**1704**  
22 cm  
8 3/4"



**1705**  
Antkoviak



**1706**  
Ermert



**1707**  
Ermert

## 04 SECTION

### Wax Knives



**1708**  
Fahnenstock  
Fig. 1  
12,5cm



**1709**  
Fahnenstock  
Fig. 2  
17cm



**1710**  
Lessmann  
Fig. 1  
12,5cm



**1711**  
Lessmann  
Fig. 2  
17cm



**1712**  
Con scodellino  
Fig. 1  
12,5cm



**1713**  
Con scodellino  
Fig. 2  
17cm



**1714**  
Fig. 1  
17 cm



**1715**  
Gritman  
Fig. 2  
17 cm




**1716**  
Fig. 3  
16 cm



**1717**  
Fig. 4  
16 cm

**Explorers**





**1808**  
Fig.1  
Spreader

**1809**  
Fig.2  
Spreader

**1810**  
Fig.3  
Spreader


**1811**  
Fig.4  
Spreader

**1812**  
Fig.5  
Spreader

**1813**  
Fig.6  
Spreader

**1814**  
Fig.7  
Spreader

### Double-Ended Root Canal Explorers Plugger/Spreaders



Root-canal pluggers

Ø 0.5    Ø 0.7    Ø 0.9    Ø 1.1



**1827**  
Fig.1  
Luks

**1828**  
Fig.2  
Luks

**1829**  
Fig.3  
Luks

**1830**  
Fig.4  
Luks



**1831**  
Fig.1  
Vertical  
condenser



**1832**  
Fig.2  
Vertical  
condenser



**1833**  
Fig.3  
Vertical  
condenser



**1834**  
Fig.4  
Vertical  
condenser


Endodontic set



**1835**  
Measuring Plate

### Root-canal pluggers

Optimum Flexibility  
Very Fine Points



Root-canal explorers



**1846**  
Fig. 6/1



**1847**  
Fig. 6/17



**1848**  
Fig. 6/23

Root-canal excavators



**1849**  
Fig. 31L



**1850**  
Fig. 32L



**1851**  
Fig. 33L

## 04 SECTION

### Excavators

1 mm



**1852**  
Fig 5

1.2 mm



**1853**  
Fig 14

1 mm



**1854**  
Fig 17

1.2 mm



**1855**  
Fig 18

1.6 mm



**1856**  
Fig 19

2 mm



**1857**  
Fig 20

1 mm



**1858**  
ASH  
Fig.243

1,4 mm



**1859**  
ASH  
Fig.G2

0,8 mm



**1860**  
Fig.1

0,8 mm



**1861**  
Fig.2

1 mm



**1862**  
Fig.11

1 mm



**1863**  
Fig.12

1,2 mm



**1864**  
Fig.21

1,2 mm



**1865**  
Fig.22

1,6 mm



**1866**  
Fig.31

1,6 mm



**1867**  
Fig.32

2 mm



**1868**  
Fig.41

2 mm



**1869**  
Fig.42


## SECTION **05**

HANDLES FOR MOUTH MIRRORS, SPATULAS,  
NEEDLE HOLDERS, SYRINGES, DENTAL TWEEZERS,  
RUBBER DAM & CLAMPS

Page 89-112

[www.arcadia14.co.uk](http://www.arcadia14.co.uk)

# Handles for Mouth Mirrors and Mouth Mirrors





1911



1912



1913



1914




1915



1916  
Parallelometer  
mouth mirror



1917  
Parallelometer  
mouth mirror  
+ handle  
code 2101



1922  
Cone - socket




1923  
Cone - socket



1924  
Cone - socket



2000



2001



2002



2003



2004




2005




2006



2007



2008



2009



2010



2011



2012




2013




2014



## Spatulas for Plaster and Alginate







**2200**  
Mayo-Hegar  
16 cm



**2201**  
Mayo-Hegar TC  
16 cm



**2202**  
Mayo-Hegar  
18 cm



**2203**  
Mayo-Hegar TC  
18 cm



**2204**  
Mini-Ryder  
18 cm



**2205**  
Mini-Ryder  
15 cm



**2206**  
Derf  
12 cm



**2207**  
Hegar-Baumgarther  
14 cm

# 05 SECTION

## Needle Holders



**2208**  
Crile-Wood  
15 cm



**2209**  
Crile-Wood TC  
15 cm



**2210**  
Adson  
17.5 cm



**2211**  
Olsen-Hegar  
14 cm



**2212**  
Baumgarther  
14.5 cm



**2213**  
Gillies  
16 cm



**2214**  
Gillies  
16 cm



**2215**  
Hegar-Baumgarther  
14 cm



**2216**  
Toennis  
18 cm



**2217**  
Toennis  
18 cm



**2218**  
Lichtenberg  
17 cm



**2219**  
Lichtenberg  
17 cm



**2220**  
Eiselberg-Mathieu  
19.5 cm



**2221**  
Boynton  
12 cm



**2222**  
Arruga  
16 cm



**2223**  
Boynton  
12 cm



**2224**  
Arruga  
16 cm

# 05 SECTION

## Needle Holders



**2225**  
Castroviejo  
14 cm



**2226**  
Castroviejo  
14 cm



**2227**  
Castroviejo  
14/60 cm



**2228**  
Castroviejo  
18 cm



**2229**  
Castroviejo  
18 cm



**2230**  
Castroviejo  
14 cm



**2231**  
Castroviejo  
14 cm



**2232**  
Castroviejo  
14/60 cm



**2233**  
Castroviejo  
18 cm



**2234**  
Castroviejo  
18 cm



**2235**  
Mathieu  
14 cm



**2236**  
Mathieu  
14 cm



**2237**  
Mathieu  
14 cm



**2238**  
Mathieu  
14 cm



**2239**  
Mathieu  
17 cm



**2240**  
Mathieu  
17 cm



**2241**  
Mathieu  
14 cm



**2242**  
Micro-Barraquer  
18 cm



The tips have tungsten carbide inserts for better handling.



**2243**  
Mathieu  
17 cm



**2244**  
Mathieu  
17 cm



**2245**  
Micro-Barraquer  
14 cm



**2246**  
Castroviejo  
13 cm



**2247**  
Castroviejo-Barraquer  
13 cm



**2248**  
Castroviejo  
13 cm



**2249**  
Castroviejo  
13 cm



**2251**  
18.03mm, cm



**2252**  
18.04mm, cm



**2253**  
18.05mm, cm




**2250**  
Fig. 6FR  
Fraizer  
15.02mm, cm



**2254**  
Fig. 12FR  
Fraizer  
15.04mm, cm



**2255**  
Fig. 10FR  
Fraizer  
15.03,3mm, cm



**2256**  
Fig. 8FR  
Fraizer  
15.02,7mm, cm



**2300**

Water Syringe complete with  
1 cannula (Luer-Lock)



**2301**

Water Syringe complete with  
1 cannula (Luer-Lock)



**2302**

Metal water syringe, self-filling,  
with 1 cannula (Luer-Lock)



**2303**

Sparrel Barrel with  
plunger suitable for  
syringes



**2304**

Cannula, only, with  
Luer-Lock connector  
suitable for syringes

## 05 SECTION

### Syringes, Rulers, Band Adaptor, Position Gauge



**2305**  
Aspirating Syringe



**2306**  
Position Gauge



**2307**  
Breech Loading Syringe



**2308**  
Aspirating & Breech Syringe



**2309**  
Band Adaptor



**2310**  
Endo Finger Ruler  
With Adjustable Slide



**2311**  
Aspirating Syringe



**2312**



**2313**



**2314**



**2315** 1,8 ml  
**2316** 2,2 ml

Cartridge Syringe complete  
With 2 tips and without needles



**2317**  
Tip Only, short



**2318**  
Tip Only, long



**2319**



**2320**



**2321**



**2322**  
Cork screw tip



**2323**  
Cork screw tip



**2324**  
Folding type, cork  
screw tip



**2325**  
Folding type, cork  
screw tip



**2326**



**2327**  
Cork screw tip



**2328**  
Spare cap with  
metric thread



**2329**  
With Whitworth  
thread



**2330**


Autoclavable Plastic Protection sleeve



**2331**  
Brass syringe with stainless steel moving parts




**2332**  
Brass syringe with stainless  
steel moving parts




smooth

**2400**  
College  
15 cm




serrated

**2401**  
College  
15 cm




smooth

**2402**  
College  
15 cm




serrated

**2403**  
College  
15 cm




smooth

**2404**  
College  
15 cm




serrated

**2405**  
College  
15 cm




smooth

**2406**  
College  
15 cm




serrated

**2407**  
College  
15 cm




smooth

**2408**  
Meriam  
16 cm




serrated

**2409**  
Meriam  
16 cm




smooth

**2410**  
Meriam  
16 cm




serrated

**2411**  
Meriam  
16 cm




smooth

**2412**  
Meriam  
16 cm



serrated

**2413**  
Meriam  
16 cm




smooth

**2414**  
Meriam  
16 cm



serrated

**2415**  
Meriam  
16 cm




**2416**  
College  
14.5 cm, 5 3/4"

**2417**  
Flagg  
16 cm, 6 1/4"

**2418**  
Flagg  
16 cm, 6 1/4"

**2419**  
Landgraf  
15 cm, 6"

**2420**  
15.5 cm, 6 1/8"



**2421**  
15.5 cm, 6 1/8"

**2422**  
London-College  
15 cm, 6"

**2423**  
London-College  
15 cm, 6"

**2424**  
London-College  
15 cm, 6"

**2425**  
London-College  
15 cm, 6"



**2426**  
CK1  
15 cm, 6"



**2427**  
CK2  
15 cm, 6"



**2428**  
GF2  
16 cm, 6 1/4"



**2429**  
Gf1  
16 cm, 6 1/4"



**2430**  
Kopp  
15.5 cm, 6 1/8"



**2431**  
Kopp  
15.5 cm, 6 1/8"



**2432**  
15 cm, 6"




**2433**  
16.5 cm, 6 1/2"



**2434**  
15.5 cm, 6 1/8"



**2435**  
15.5 cm, 6 1/8"



**2436**  
London-College  
15 cm, 6"

**2437**  
London-College  
15 cm, 6"

**2438**  
Meriam  
16 cm, 6 1/4"

**2439**  
Allen  
18 cm, 7 1/8"




**2440**  
Albrecht  
16 cm, 6 1/4"

**2441**  
Albrecht  
16 cm, 6 1/4"


**2442**  
16 cm, 6 1/4"

**2443**  
16 cm, 6 1/4"

**2444**  
16 cm, 6 1/4"



**CLAMPS ( PLIERS )**



# Rubber Dam Clamps

**Anterior**



**2510 - 9**  
Winged  
Labial Clamp



**2511 - 9T**  
Winged  
Edged Labial Clamp



**2512 - 212**  
Winged  
Cervical Clamp



**2513 - 210**  
Winged  
Labial Clamp

**Premolars**



**2514 - 00**  
Winged



**2515 - 0**  
Winged



**2516 - W0**  
Wingless

For Small Teeth, Elongated or With Narrow Collars



**2517 - 1**  
Winged



**2518 - 1T**  
Winged



**2519 - W1**  
Wingless

For General Use With Adjustable Jaws According  
To Gingival Line Shapes



**2520 - 2A**  
Winged



**2521 - 2AT**  
Edged, Winged



**2522 - W2A**  
Wingless

For Premolars, With Elongated Collar



**2523 - 2**  
Winged



**2524 - 2T**  
Edged, Winged



**2525 - W2**  
Wingless

For Premolars, With Elongated Collar



**2526 - 29**  
Wingless



**2527 - 60**  
Wingless



**2528 - 3**  
Winged



**2529 - W3**  
Wingless

**For Small Lower Molars**



**2530 - 4**  
Winged



**2531 - W4**  
Wingless

**For Small Lower Molars With Round clamps**



**2532 - 7**  
Winged



**2533 - W7**  
Wingless

**For Square Lower Molars**



**2534 - 8**  
Winged



**2535 - W8**  
Wingless

**For Round Upper Molars**



**2536 - 8A**  
Winged



**2537 - W8A**  
Wingless

**For Molar Roots**



**2538 - 12A**  
Winged



**2539 - 13A**  
Wingless

**For Partially Broken Molars**



**2540 - 14**  
Winged



**2541 - 14T**  
Edged Winged



**2542 - W14**  
Wingless



**2543 - 14A**  
Winged



**2544 - W14A**  
Wingless

**For Round Molars, Partially Broken, Undersized or Inclined**



**2545 - 26**  
Wingless



**2246 - 51**  
Wingless



**2547 - 202**  
Winged



**2548 - 205**  
Winged

## SECTION **06**

MATRIX RETAINER, MATRIX BANDS,  
CROWN REMOVERS, SCISSORS, FORCEPS,  
AMALGUM GUNS, CARRIES, COPPER RING REMOVING  
PLIERS, PAPER ARTICULATOR FORCEPS

Page 113-124

[www.arcadia14.co.uk](http://www.arcadia14.co.uk)





**2617**  
**Molar**  
**Fig 7**



**2618**  
**Molar**  
**Fig 8**



**2619**  
**Molar**  
**Fig 9**



**2620**  
**Molar**  
**Fig 10**



**2621**  
**Tofflemire**



**2622**  
**Tofflemire**



**2623**  
**Tofflemire**



**2624**  
**Bicuspid**



**2625**  
**Bicuspid**



**2626**  
**Bicuspid**



**2627**  
**Bicuspid**



**2628**  
**Bicuspid**



**2629**  
**Bicuspid**



**2633**




**2630**  
**Molar**




**2631**  
**Molar**



**2632**  
**Molars**



**2704**  
Crown Spreader, Crown Remove.  
Special instrument for crown spreading  
after separating, simple and safe  
operation.





Universal "HM"  
12 cm, 4 3/4"



2709



2710



2711  
Beebe  
10,5 cm, 4 1/8"



2712  
Merston



2713  
Morell



2714  
Morell



2715  
Morell



2716  
Morell



2717  
Treyman



2718  
Automatic Crown Remover  
With 2 Points



2719



2720



**2810**  
Amalgam Well  
Non Skid Base  
Ø 3.5 mm

**Copper Ring Removing Pliers, Paper Articulator Forceps**

**Copper ring removing pliers**



**2900**  
Bohm  
11,5 cm



**2901**  
Bohm  
14,5 cm



**2902**  
Furrer  
16 cm



**2903**  
Furrer  
16 cm



**2904**  
Miller  
15 cm



**2905**  
Miller  
15 cm



**2906**  
Kelly  
14 cm



**2907**  
Halstead-Mosquito  
12 cm



**2908**  
Micro-Mosquito  
10 cm



**3000**  
Friedmann  
14 cm, 5 1/2"



**3001**  
14 cm, 5 1/2"



**3002**  
Mead  
14 cm, 5 1/2"



**3003**  
Mead  
14 cm, 5 1/2"



**3004**  
Mead  
17 cm, 6 3/4"



**3005**  
Schell  
17 cm, 6 3/4"



**3006**  
Luer  
15 cm



**3007**  
Friedman  
14 cm



**3008**  
Mini Friedman  
12 cm



**3009**  
Goldman-Fax  
10 cm



**3010**  
Blumenthal  
15 cm



**3011**  
Blumenthal 45  
14,5 cm



**3012**  
Blumenthal 90  
14 cm



**3013**  
Bohler  
15 cm



**3014**  
Kerrison  
15 cm



**3015**  
Blumenthal  
16 cm, 6 1/4"



**3016**  
Blumenthal  
16 cm, 6 1/4"



**3017**  
Blumenthal  
16 cm, 6 1/4"



**3018**  
Kleinert-Kutz  
14 cm, 5 1/2"



**3019**  
Ruskin-Mini  
15 cm



**3020**  
Beyer  
18 cm, 7 1/8"



**3021**  
Zaufal-Jansen  
18 cm, 7 1/8"



**3022**  
Niro  
No. 5 S, 14 cm, 5 1/2"



**3023**  
Cleveland  
No. 5 S, 17 cm, 6 3/4"



**3024**  
No. 18, 15 cm, 6"



**3025**  
Cleveland  
16 cm, 6 1/4"




**3026**  
Gardner  
16 cm, 6 1/4"



**3027**  
Liston-mini  
15 cm, 6"


# 06 SECTION

## Bone Elevators



**3100**  
Williger  
16 cm, 6 1/4"

**3101**  
Henahan  
17 cm, 6 3/4"




**3102**  
Freer  
19 cm, 7 1/2"

**3103**  
Bennet  
17 cm, 6 3/4"

**3104**  
Seldin  
19 cm, 7 1/2"

**3105**  
Seldin  
19,5 cm, 7 3/4"



## SECTION **07**

**PLIERS FOR ORTHODONTICS  
& IMPRESSION TRAYS**

**Page 125-139**

**[www.arcadia14.co.uk](http://www.arcadia14.co.uk)**

Pinze per formare loops  
Loops forming pliers



**3200**  
Adams  
Adams pliers for bending  
wires up to 0,7mm. (.28").



**3201**  
Adams  
Adams pliers for bending  
wires up to 0,7mm. (.28").



**3202**  
Adams  
Tweed loop forming pliers.  
For forming loop in wires up  
to 0,56 mm. (.22").



**3203**  
Adams  
Aderer three jaws bending  
pliers. For wires up  
to 0,7 mm. (.28").



**3204**  
Adams  
Aderer three jaws bending  
pliers. For wires upto 0,7 mm.  
(.28").

Ligature Cutter Pliers



**3205**

Micro-miniature light wire cutter. It cuts soft ligature wires. For use in caso of limited access. Caution: maximum cutting capacity: soft ligature wire up to 0,30 mm. (.012").



**3206**

Wire cutter. For ligature wires and pins. It is strong enough to bend archwire ends after distal cutting. Maximum cutting capacity up to 0,56 x 0,7 mm. (.022" x .028").



**3207**

Wire cutter. For ligature wires and pins. It is strong enough to bend archwire ends after distal cutting. Maximum cutting capacity: up to 0,9 mm. (.036").



**3208**

Wire cutter. For ligature wires and pins. It is strong enough to bend archwire ends after distal cutting. Maximum cutting capacity: up to 1,5 mm. (.044").

Pinze per formare laboratorio  
Laboratory pliers



**3209**

Flat nose pliers with smooth beaks for bending wires up to 0,7 mm. (.028").



**3210**

Bending wires up to 0,7 mm. (.028").



**3211**

Bending wires up to 0,7 mm. (.028").



**3212**

Bending wires up to 0,7 mm. (.028").



**3213**  
Nance

Nance loop closing pliers, for bending of stop loops and special for Adams' clasps. Wires up to 0,7 mm. (.028").



**3214**  
Nance

Universal pliers for bending wires up to 0,9 mm. (.035") or cutting up to 0,7 mm. (.028").



**3215**

Stainless steel pliers tray  
standard dimensions



**3216**  
Universal

Universal pin roughening, cutting and bending pliers. For wires up to 0,9 mm. (.035"). Will cut wires up to 0,7 mm. (.028").

**Utility and specialty pliers**



**3217**  
Reynolds

Reynolds contouring pliers.



**3218**  
Reynolds

Lingual arch pliers used to form double and triple back bends in 0,7 mm. (.028") or 0,9 mm. (.035") wire for insertion in lingual sheaths.



**3219**  
Weingart

Weingart pliers with serrated points.



**3220**  
De La Rosa

De La Rosa arch forming pliers. Designed to form and contour arch wires on round or rectangular wire. The beak has grooves of 0,4 mm (.016"), 0,46 mm (.018") and 0,56 mm (.022").



**3221**  
De La Rosa

Arch forming pliers with a non grooved contouring surface the forming of precise arches.

### Loop formare pliers



**3222**  
Waldsachs

Waldsachs universal pliers for bending wires up to 1,0 mm. (.039") or cutting up to 0,7 mm. (.028").



**3223**

Round and Concave pliers for loop forming and bending of thin round and square wires up to 0,5 mm. (.020").



**3224**  
Tweed

Tweed loop Forming Pliers. For round or rectangular wires up to 0,56 x 0,7 mm. (.022" x .028"). Replaceable cylindrical beak with three sizes: 1,15 mm. (.045), 1,50 mm. (.059"), 1,90 mm. (.075"). Concave beak with slight parallel serrations, to keep the wire in proper position.



**3225**  
Tweed

Tweed Loop Forming Pliers. For making perfect Omega Loops with round or rectangular wires up to 0,56 x 0,7 mm. (.022" x .028"). The concave beak is exactly the same as pliers 5809-1



**3226**  
Tweed

Single diameter loop forming pliers. Wire up to 0,56 x 0,7 mm. (.022" x 0.28"). Replaceable concave beak has only one diameter section of 1,50 mm. (.059").

**Distal end cutters**



**3227**

Distal end cutter -safety hold- a precision cutter designed for distal arch cutting in the patient's mouth. It cuts hard wires and holds the loose distal end preventing any damaging of oral tissue. Maximum cutting capacity rectangular wires up to 0,56 x 0,7 mm. (.022" x .020").



**3228**

Micro- distal end cutter - safety hold-. Similar to 3227



**3229**

Direct bonding bracket remover (wide). The removal grip firmly for easy incisal-gingival (vertical) use or, if needed, for mesial-distal (horizontal) use with the same ease.



**3230**

Direct bond attachment remover (angled). For removing all bonded attachments. May be used mesial-distal-horizontal) or incisal-gingival (vertical) in anterior areas. Angled perfectly for easy access to difficult posterior areas.



**3231**

Wire bending pliers. For wires up to 0,9 mm. (.036").



**3232**

Bird Beak Angle pliers for bending round, square or rectangular wire up to 0,7 mm. (.028").



**3233**

Double rounded jaw wire bending pliers. For wires up to 0,7 mm. (.028").



**3234**

Bird beak pliers for bending round, square or rectangular wire up to 0,7 mm. (.028"). With cutter.






**3235**

Double rounded jaw wire bending pliers. For wires up to 0,7 mm. (.028").

**Bird beack pliers**

# Pliers for Othodontics and Prosthetics





**3251**  
Langenbeck  
140 mm



**3252**  
Langenbeck  
130 mm



**3253**  
Langenbeck  
130 mm



**3254**  
Peeso  
125 mm



**3255**  
Goslee  
125 mm



**3256**  
Fischer  
140 mm



**3257**  
Andresen  
130 mm



**3258**  
Schwarz  
130 mm



**3259**  
MSY  
130 mm



**3260**  
Langenbeck  
125 mm



**3262**



**3263**



**3264**



**3265**



**3266**



**3267**



**3268**

# Pliers for Othodontics and Prosthetics




3269



3270



3271



3272



3273



3274



3275



3276



3277



3278



3279



3280



3281



3282



3283



3284



3285



3286

All Impression Trays Upper & Lower Type  
are available in S, M, L & XL Sizes



**3300**  
"HYDRO-LOCK"  
stainless steel  
water-cooled  
impression trays  
with retention rim.



**3301**  
"HYDRO-LOCK"  
stainless steel  
water-cooled impression  
trays with retention rim.



**3302**  
"NEW SUPER"  
stainless steel  
impression trays  
with retention rim.



**3303**  
"NEW SUPER" stainless  
steel impression trays  
with retention rim.



**3304**  
"ANATOMIC" solid  
regular stainless  
steel impression  
trays without  
retention rim.



**3305**  
"ANATOMIC" solid regular  
stainless steel impression  
trays without retention rim.



**3306**  
Perforated-total denture  
stainless steel  
impression trays with rim.



**3307**  
Perforated-total denture  
stainless steel  
impression trays with rim.

## Stainless Steel Impression Trays

All Impression Trays Upper & Lower Type  
are available in S, M, L & XL Sizes



**3308**  
"PERMA-LOCK" solid regular stainless  
steel impression trays with retention rim.



**3309**  
"PERMA-LOCK" solid regular stainless  
steel impression trays with retention rim.



**3310**



**3311**



**3312**



**3313**



**3314**



**3315**

### PACKINGS



**3316**  
Set of six



**3317**  
Set of Eight



**3318**  
Solid stainless steel  
anterior depressed  
impression trays with  
retention rim.



**3319**



**3320**



**3321**



**3322**



**3323**



**3324**



**3325**  
Solid stainless steel  
anterior depressed  
impression trays with  
retention rim.

**All Impression Trays Upper & Lower Type  
are available in S, M, L & XL Sizes**



**3326**  
Fig. 0 per bambini - for children



**3327**  
For narrow toothed upper jaws



**3328**  
Fig. 0 per bambini - for children



**3329**  
For partially toothed upper jaws



**3330**  
For partially toothed lower jaws



**3331**  
For edentulous upper jaws



**3332**  
For edentulous lower jaws

## Stainless Steel Impression Trays

All Impression Trays Upper & Lower Type  
are available in S, M, L & XL Sizes



**3333**  
Solid total denture stainless steel impression  
tray with rim.



**3334**  
Solid total denture stainless steel impression  
tray with rim.



**3335**  
"PERMA-LOCK" perforated regular stainless  
steel impression trays with retention rim.



**3336**  
"PERMA-LOCK" perforated regular stainless  
steel impression trays with retention rim.



**3337**  
Perforated stainless steel anterior depressed  
impression trays with retention rim.



**3338**



**3339**



**3340**



**3341**



**3342**



**3343**



**3344**  
Perforated stainless steel anterior depressed  
impression trays with retention rim.



**3345**  
Perforated stainless  
steel impression trays  
without retention rim  
for complete denture.



**3346**  
Perforated stainless steel  
impression  
trays without retention  
rim for complete denture.

## SECTION **08**

ARRICULATORS, DENTAL  
& SURGICAL HOLLOWARES

Page 140-144

[www.arcadia14.co.uk](http://www.arcadia14.co.uk)



3400



3401



3402



3403



3404



3405



3406



3407



3408



3409



**3500**  
**ENDODONTIC BOX**



**3501**  
**STERILIZING BOX**  
Available in dif. sizes  
17x17x3 cm  
20x10x4 cm  
25x12x6 cm  
30x14x6 cm



**3502**  
**ENDO TRAY BOX**



**3503**  
**COTTON DISPENSER**



**3504**  
**ALCOHOL LAMP**

**Sterilizing Drums**



**3505**



**3506**



**3507**

**Sterilizing Boxes**



3508



3509



3510

**Sterilizing Boxes**



3511



3512



3513



3514



3515



3516

**FORCEP AND THERMOMETER JARS**

**COTTON DISPENSER**



3517  
IRRIGATOR(Douch can)



3518  
MEASURING JUG



3519  
KIDNEY TRAY



**3520**  
**FORCEP JAR**



**3521**  
**GALLIPOT WITH & WITHOUT LID**



**3522**  
**SPITTEN MUG**



**3523**  
**SPIRIT LAMP**



**3524**  
**Perforated Tray**



**3525**  
**Perforated Base**



**3526**  
**Solid Base**



**3527**  
**Instruments Tray**



**3528**  
**Instruments Tray**



**3529**  
**Instruments Tray  
(Perforated)**



**3530**  
**Instruments Tray Shallow  
(Mayo)**



**3531**  
**Scalers Tray**



**3532**  
**Cathetar Tray**

# Alphabetical Index

|                        |             |                     |                     |                           |             |
|------------------------|-------------|---------------------|---------------------|---------------------------|-------------|
| <b>A</b>               |             |                     | <b>G</b>            |                           |             |
| Adams                  | 125, 133    | Cleveland           | 122                 | Goslee                    | 134         |
| Aderer                 | 133         | College             | 105, 107            | Gouges                    | 68          |
| Adson                  | 96          | Columbia            | 64, 66, 74          | Gracey Curettes           | 69          |
| Ainsworth              | 110         | Columbia            | 79                  | Gracey Mesial Distal      | 69          |
| Albee                  | 77          | Cone-Socket         | 90                  | Gracey                    | 56-69-72    |
| Albrecht               | 109         | Cork Screw Tip      | 104                 | Gritman                   | 76-78-80    |
| Alcohol Lamp           | 142         | Cotton Dispenser    | 142, 143            |                           |             |
| Alginate               | 93-94       | Coupland            | 37                  | <b>H</b>                  |             |
| Allen                  | 109         | Crane-kaplan        | 74                  | Halstead-Mosquito         | 118         |
| Amalgam Carriers       | 117         | Crile-Wood          | 96                  | Hegar-Baumgartner         | 95-96       |
| Amalgam Guns           | 117         | Crown Removers      | 115-116             | Heidbrink                 | 45-48-50    |
| Amalgam Well           | 117         | Cryer               | 35, 37, 48          | Heidemann                 | 55          |
| Andresen               | 134         | Curette             | 64, 66              | Hemingway                 | 75          |
| Angle                  | 133         |                     |                     | Henahan                   | 123         |
| Antkoviak              | 78          | <b>D</b>            |                     | How                       | 133         |
| Apical                 | 36, 44      | De La Rosa          | 129                 | Hylin                     | 40, 76      |
| Arruga                 | 97          | Dental Probes       | 57                  |                           |             |
| Articulators           | 140         | Dental Syringes     | 103-104             | <b>I</b>                  |             |
| Ash                    | 53, 88      | Derf                | 95                  | Impression Trays          | 136-139     |
| Aspirating Syringes    | 102         | Didner              | 113                 | Intra-Ligamental Syringes | 104         |
|                        |             |                     |                     | Irrigator                 | 143         |
| <b>B</b>               |             | <b>E</b>            |                     | Ivory                     | 110, 113    |
| Band Adaptor           | 102         | Eliselberg-Mathieu  | 97                  | Iwanson                   | 77          |
| Barry                  | 47, 50      | Elliott             | 53                  |                           |             |
| Baumgarthner           | 96          | Endo Finger Ruler   | 102                 | <b>J</b>                  |             |
| Beale                  | 76          | Endo Tray Box       | 142                 | Jaquette                  | 74          |
| Beebee                 | 116         | Endodontic Box      | 142                 | Johnson                   | 133         |
| Bein                   | 33, 40      | Ermert              | 78                  |                           |             |
| Bennet                 | 53, 123     | Excavators          | 87-88               | <b>K</b>                  |             |
| Beyer                  | 121         | Explorers           | 81                  | Kelly                     | 118         |
| Bicuspids              | 114         | Extracting Forceps  | 01-31               | Kerrison                  | 120         |
| Birdbeak Pliers        | 132         |                     |                     | Kidney Tray               | 143         |
| Blumenthal             | 120, 121    | <b>F</b>            |                     | Kilian                    | 124         |
| Bohler                 | 120         | Fahnenstock         | 79                  | Kilian-Claus              | 124         |
| Bohm                   | 118         | Filling Instruments | 52-56               | Klein                     | 07          |
| Boley                  | 77          | Fischer             | 134                 | Kleinert-Kutz             | 121         |
| Bone Chisels           | 68          | Flagg               | 107                 | Kopp                      | 39, 49, 108 |
| Boynton                | 97          | Flohr               | 38                  |                           |             |
| Breech Loading Syringe | 102         | Forceps Jar         | 143, 144            | <b>L</b>                  |             |
| Brewer                 | 110         | Fox                 | 58                  | Landgraf                  | 107         |
| Buck                   | 61-63, 65   | Fraizer             | 100                 | Langenbeck                | 134         |
| Buckley                | 67, 124     | Freer               | 123                 | Le Cron                   | 76          |
| Burnishers             | 58          | Friedman            | 119-120             | Lecluse                   | 43          |
|                        |             | Furrer              | 119                 | Lessmann                  | 79          |
| <b>C</b>               |             |                     |                     | Lichtenberg               | 97          |
| Caroll                 | 77          | <b>G</b>            |                     | Listion-mini              | 122         |
| Castroviejo            | 77, 98, 100 | Gallipot            | 144                 | London Hospital           | 41          |
| Cathetar Tray          | 144         | Gardner             | 67, 122             | London-College            | 107, 109    |
| Cement Spatulas        | 91, 92      | GF1                 | 108                 | Lucas                     | 75          |
| Chompret               | 44          | GF2                 | 108                 | Luer                      | 120         |
| Ck1                    | 108         | Gillies             | 96                  | Luks                      | 84          |
| CK2                    | 108         | Goldman Fox         | 60, 63, 65, 73, 120 |                           |             |

## M

|                      |             |
|----------------------|-------------|
| Machtou              | 83          |
| Mallets              | 68          |
| Mathieu              | 99, 100     |
| Matrix Bands         | 113, 114    |
| Matrix Retainers     | 113         |
| Mayo-Hegar           | 95          |
| Mc Call              | 73          |
| Mc Kellops           | 133         |
| Mead                 | 43, 58, 119 |
| Measuring Jug        | 143         |
| Measuring Plate      | 84          |
| Medicon              | 78          |
| Meriam               | 106, 109    |
| Merrit A             | 60          |
| Mershon              | 110, 116    |
| Micro-Barraquer      | 99-100      |
| Micro-Mosquito       | 118         |
| Miller               | 118, 124    |
| Mini-Ryder           | 95          |
| Molar                | 114         |
| Morell               | 116         |
| Mouth Mirror Handles | 89, 90      |
| Mouth Mirrors        | 89, 90      |
| MSY                  | 76, 134     |

## N

|        |          |
|--------|----------|
| Nabers | 83       |
| Nance  | 128, 133 |
| Niro   | 122      |

## O

|             |            |
|-------------|------------|
| Obwegeser   | 124        |
| Ochsenbein  | 61, 62, 67 |
| Olsen-Heger | 96         |
| Oms         | 83         |
| Ott         | 41         |

## P

|                  |       |
|------------------|-------|
| P                | 76    |
| P.K. Thomas      | 68    |
| Partsch          | 134   |
| Pocket Probes    | 58-60 |
| Plaster Knives   | 80    |
| Plaster Spatulas | 93-94 |
| Position Gauge   | 102   |
| Pott             | 42    |

## R

|          |        |
|----------|--------|
| Rein     | 81     |
| Relax    | 05     |
| Reynolds | 129    |
| Rhodes   | 61, 62 |

|                       |         |
|-----------------------|---------|
| Rhodium               | 89      |
| Roach                 | 76      |
| Root Elevators        | 32-51   |
| Root Splinter         | 31      |
| Root-Canal Excavators | 86      |
| Root-Canal Explorers  | 86      |
| Root-Canal Plungers   | 84-85   |
| Rubber Dam Clamps     | 111-112 |
| Rubber Dam Frame      | 110     |
| Ruskin-Mini           | 121     |

## S

|                   |             |
|-------------------|-------------|
| Scalers           | 58          |
| Schell            | 119         |
| Schilder          | 85          |
| Schluger          | 61, 62      |
| Schwarz           | 134         |
| Seldin            | 34, 48, 123 |
| Spirit Lamp       | 144         |
| Spitten Mug       | 144         |
| Spreader          | 82-83       |
| Sterilizing Boxes | 142, 143    |
| Sterilizing Drums | 142         |
| Sugarman          | 61, 162     |

## T

|                      |     |
|----------------------|-----|
| Taylor               | 74  |
| Toennis              | 97  |
| Tofflemire           | 114 |
| Tofflemire-Junior    | 113 |
| Tofflemire-Senior    | 113 |
| Tofflemire-Universal | 113 |
| Towner Jaquette      | 74  |
| Towner               | 74  |
| Trays                | 144 |
| Treymann             | 116 |
| Tweed                | 130 |

## U

|                        |          |
|------------------------|----------|
| Univarsal "HM"         | 116      |
| Univarsal              | 128, 133 |
| University of Michigan | 60       |
| USA Model              | 133      |

## V

|                    |    |
|--------------------|----|
| Vehe               | 76 |
| Vernier            | 77 |
| Vertical Condenser | 83 |

## W

|               |            |
|---------------|------------|
| Wakefield     | 67         |
| Waldron       | 133        |
| Waldsachs     | 130, 133   |
| Warwick-James | 41         |
| Water Syringe | 101        |
| Weingart      | 129, 133   |
| White         | 52, 60, 83 |
| Williams      | 123        |
| Williger      | 46, 50     |
| Winter        | 54, 55     |
| Woodson       |            |

## Y

|              |     |
|--------------|-----|
| Young        | 133 |
| Younger-Good | 73  |

## Z

|                |        |
|----------------|--------|
| Zahle          | 76     |
| Zaufal-Jansen  | 121    |
| Zerfing-Chisel | 64, 66 |
| Zuricher       | 77     |



**MANUFACTURER**



**SEHAR BATOOL INTERNATIONAL**

S.I.E Fateh Garh Near Umer Town Sialkot, Pk

Tel: +92 52-3552898 Fax: +92 52-3251597

Email: [seharbatool@brain.net.pk](mailto:seharbatool@brain.net.pk) / [Info@seharbatool.com](mailto:Info@seharbatool.com)

Website: [www.seharbatool.com](http://www.seharbatool.com)